

ADDING A NEW DRUG INTO CONSIS

1	CREATE STOCKCARD IN THE CONSIS SOFTWARE	→ Select CREATE PRODUCT
		→ Select the PRODUCT CODE field
		→ Scan the barcode of the new product to be added
		→ Ensure the following fields are complete for the new drug: Product Name, Product length (mm), Dosage Form, Package Unit
		*Tick the EXCESS STOCK field if the item will typically have excess stored in the pharmacy e.g., those items ordered in bulk quantities
		HINT - You may wish to proceed directly to step 3 'Add channel allocation to new drug' by selecting ADD CHANNEL or select OK to save
2	PREPARE AND CHOOSE AN APPROPRIATE CHANNEL	Choose an appropriate channel from the list of free channels (ensuring correct package alignment). The free channel list may be filtered according to height by ticking or unticking the channel colour boxes; Red, Green, Blue or Yellow. → Adjust dividers / bottle channels to ensure the new drug fits snugly and is centred over the pick gap.
3	ADD CHANNEL ALLOCATION TO NEW DRUG	→ Select CREATE PRODUCT
		→ Look up drug by either scanning (select PRODUCT CODE field) or you may prefer typing the name, if so (select PRODUCT SEARCH field)
		→ Select EDIT to enter the product stock card
		→ Select ADD CHANNEL
		→ Select OK to save
5	FILL DRUG INTO CHANNEL	→ Select STOCK PACKAGES
		→ Select DELIVERY or EXCESS
		→ Select OK
		→ Scan product
		→ Select Channel
		→ Enter quantity filled into channel → Select OK
6	ADD NEW DRUG TO DISPENSE SOFTWARE ROBOTIC LIST	NOTE - Various Dispense programs associate robotic drugs in different ways, refer to your dispense software user manual or contact their help desk. HINT - Remember to also remove the old drug from your dispense software if it is no longer in the CONSIS (check for multiple channels)
7	TEST DELIVERY THROUGH DISPENSE SOFTWARE	

If you are having trouble with your CONSIS PC
please contact the Service line
1800 94 55 22 | 8am - 8pm EST